

**COMPETENCY ASSESSMENT AND IMPROVEMENT**

Aspect of Care/Skill **PMV Observation**

Employee's Name (PRINT):

Facility: **Gaylord**

Competency Performance Check List	Able To Perform	Comments
<b>Preliminary Steps</b>		
1. Explains procedure to patient		
2. Assesses patient's ability to tolerate cuff deflation (slow, full cuff deflation and assess Vt and Phonation attempt)		
3. Proceeds with valve use if patient succeeds with step 2.		
4. Has Oral and tracheal suction available		
5. Suction Mouth and Trach		
6. Places valve in-line, slowly deflates cuff (should note PIP prior to deflation)		
7. Suction oral and trach as needed.		
8. Changes vent to NIV and turns PEEP off.		
9. Observes vent parameters-makes changes to maintain adequate Vt and PIP, increase Vt in small increments to maintain previous PIP.		
10. If in PS or SIMV set TI SPONT to 1 second.		
11. If in PC usually need to decrease PC in small increments for valve tolerance.		
12. In NIV 840 allows for low pressure alarm-set to at least 10ccmH2o.		

DATE: \_\_\_\_\_ CHECK ONE:  PASSED  NEEDS TO REPEAT

COMPETENCY VALIDATED BY (observer):

PRINT NAME: \_\_\_\_\_

SIGNATURE: \_\_\_\_\_